



## SINGLE and SPU

---



**EASPU with 3PS**  
Easy Applicator Single Patient Use

## EAS - Easy Applicator Single - Product Overview



The **EAS** is as Easy to Apply and Use as a Sliding Sheet but with all the benefits of a Patient Lateral Transfer Air Mattress

**Fast and Easy to Apply & Remove**



Designed especially for use in:

**The ER - Emergency Room  
For Trauma  
For Imaging and Radiology  
For Mortuaries  
For Patients that are difficult to Log Roll**



**Can be applied without the need to Log-Roll a Patient**

### Product Highlights

- Easy to Apply and Remove - less risk of skin trauma on application and removal
- Easy to Wipe and Clean
- Patient Lateral Transfer Weight of Approximately 5% - if used as directed and instructed
- SWL (Safe Working Load) 500 Kg or 1100 Lbs
- Compatible with X-Ray and MRI Imaging
- Surface Neutral
- Waterproof
- Anti Static
- Anti Stain
- Latex Free
- High Clearance

The **EAS** has been designed as a Single Patient Use Product.

The construction and the use of specially designed and manufactured Technical Fabrics provides Superior Patient Handling and Support.

For the Care Providers, the **EAS** offers a number of benefits including easy application and removal. The **EAS** is as easy to apply and use as a Sliding Sheet (Low Friction Gliding Sheets) but with all the benefits of a Patient Lateral Transfer Air Mattress.

The **EAS** can be applied under the Patient without the need of Log-Rolling the Patient. Several techniques are available for this application process including folding and unrolling and slip/slide/pull. Your **EAS** Representative will be able to provide information and training on the various techniques.

## **EASPU - Easy Applicator Single Patient Use - Product Overview**



The **3PS System** makes Application and Transfers Easy.

The **3PS System** facilitates the Easy Application and Positioning of the Patient and helps with the Easy Placement of the Air Mattress on transfer

**With 3PS  
(Patient Positioning &  
Placement System)**



Designed especially for use in:

**The OR - Operating Room  
For Intensive and Critical Care  
For Wards**

**And as a General Purpose Patient  
Lateral Transfer Air Mattress in the  
Hospital**



### **Product Highlights**

- With a Non Slip Top Cover
- With 3PS - Patient Positioning and Placement System
- Easy to Wipe and Clean
- Patient Lateral Transfer Weight of Approximately 5% - if used as directed and instructed
- SWL (Safe Working Load) 500 Kg or 1100 Lbs
- Compatible with X-Ray and MRI Imaging
- Surface Neutral
- Waterproof
- Anti Static
- Anti Stain
- Latex Free
- High Clearance

The **EASPU** has been designed as a Single Patient Use Product.

The construction and the use of specially designed and manufactured Technical Fabrics provides Superior Patient Handling and Support.

For the Care Providers, the **EASPU** offers a number of benefits including the proprietary **3PS** - Patient Positioning and Placement System. The **3PS System** facilitates the Easy Application and Positioning of the **EASPU** and provides Easy Placement of the **EASPU** on transfer. For details and training contact your local Easy Applicator Representative.

## STANDARD STOCK SIZES OF THE SINGLE and SPU RANGE



### EAS Easy Applicator Single

SOLD IN CARTONS OF 10 UNITS



36 inches  
(Standard)



### EASPU with 3PS Easy Applicator Single Patient Use



34 inches  
(Standard)

Individually Identifiable.

Individually Traceable

The **uSN** (Product **Unique Serial Number**) on the Air Mattress provides all the information necessary about each individual product including:

- The Date of Manufacture
- Details and QC (Quality Control) of the Raw Materials utilised in the construction and manufacture of the Air Mattress
- The Operatives involved in the construction and manufacture of the Air Mattress
- The QC (Quality Control) Inspector, and
- Details of the Packing Operative.



39 inches  
(Standard Plus)



50 inches  
(Bariatric)